



P R O F E S S I O N A L

LINE FILLERS

HYALURONIC ACID

RX.HA. FILLERS
CE 0483

SWISS MADE
Be Ceuticals certified ISO 13485



RX.HA.12mg

HYALURONIC ACID

rejuvenation

12 mg/ml

RX.HA.16mg

HYALURONIC ACID

first lines



16 mg/ml

RX.HA. HYALURONIC ACID PROFESSIONAL LINE FILLERS



COMPOSITION

Monophasic densified Sodium Hyaluronate. Ideal concentrations and optimal spreading for dermal rehydration, redensification of volumes, first lines and deep folds.

INDICATIONS

RX.HA. range is a viscoelastic solution to correct fine lines or deeper wrinkles and folds, increase or restore volume to the face and to remodel facial contours. It is indicated to be injected into the deep dermis or subcutis.

Lion wrinkle correction

Orbital improvement

Cheek augmentation

Nasolabial area correction

Lips drawing

Chin augmentation


RX.HA.21mg
HYALURONIC ACID

balance



RX.HA.25mg
HYALURONIC ACID

ultra deep



GIVE YOUR PATIENT A LIFT WITH RX.HA. HYALURONIC ACID



ORBITAL IMPROVEMENT

RX.HA.16mg




LIPS DRAWING

RX.HA.21mg





NASOLABIAL AREA CORRECTION

RX.HA.25mg & 25mg+




FACIAL VOLUMES – CHIN AUGMENTATION

RX.HA.25mg+



Innovative cross linked monophasic densified sodium hyaluronate



SEE THE DIFFERENCE

RX.HA.25mg⁺

HYALURONIC ACID

intense volume



24 mg/ml


FILLERS FOR SMART DOCTORS

MONOPHASIC DENSIFIED HYALURONIC ACID

Thanks to a special process of cross-linking and purification, Be Ceuticals offers the best quality of Hyaluronic Acid gels for a unique integration into the skin and long lasting results up to 12 months.


BE CEUTICALS EXCLUSIVELY USES:

- Pure and natural HA
- 100% non-animal
- Cross-linked Monophasic Densified gels (uniforms and without particles for a better integration into the skin)
- Biofermentative origin (streptococcus equi bacteria)




THE PROCESS STEP BY STEP


A) Cross-linking




1) Chains are carefully aligned to each other.



2) Cross-linking step using BDDE Stable, very dense, three-dimensional network of HA is obtained.




3) Viscoelastic optimisation. A smooth, highly viscous and highly elastic gel is obtained for better aesthetic effect and long-lasting results.




4) Homogenization. The gel is smooth and homogenous. Injection profile is easy and regular.

C) Purification




5) All undesired compounds are removed from the gel.



6) Removal of oxygen enhances the stability of the gel once filled into the syringes.

D) Result

A highly dense three-dimensional matrix is obtained. A perfectly pure gel made of very long and strong chains of Hyaluronic Acid.



P R O F E S S I O N A L

LINE FILLERS

HYALURONIC ACID

	Hyaluronic Acid concentration	Pack size	Easyflow system	Recommended indications*
RX.HA.25mg [⊕] HYALURONIC ACID	24 mg +	2 x 2 ml	All the range: 4 x 27G	Treatment of facial volumes, intense skin depressions and deep folds
RX.HA.25mg HYALURONIC ACID	24 mg	2 x 1 ml		Treatment of deep skin depressions and cheek augmentations
RX.HA.21mg HYALURONIC ACID	20 mg	2 x 1 ml		Treatment of moderate skin depressions and lips drawing
RX.HA.16mg HYALURONIC ACID	16 mg	2 x 1 ml		Treatment of light skin depression and orbital area
RX.HA.12mg HYALURONIC ACID	12 mg	2 x 2 ml		Hydration and rejuvenation. Non-crosslinked Hyaluronic Acid

YOUR SOLUTION
FROM FIRST LINES
TO DEEP INTENSE SKIN
DEPRESSIONS



*For more information, see individual product package insert



SWISS MADE
Be Ceuticals certified ISO 13485

FILLERS FOR SMART DOCTORS



┌
RX.HA.
25 mg **+**
2ml FILLERS
└



BE CEUTICALS SA
CH-1964 CONTHEY - SWITZERLAND
be-ceuticals.swiss

SWISS MADE
Be Ceuticals certified **ISO 13485**

